



endoit[®]

WIRELESS ENDODONTIC HANDPIECE



Torque range : 0.5 ~ 4.0 N.cm

It can be adjusted from 0.5 Ncm to maximum 4.0 N.cm and can be set in 0.5 N.cm units.



RPM range : 120 ~ 500rpm

The rotation speed can be set from 120 rpm to 500 rpm.



5 Operation Modes

Forward, Reverse, Auto Reverse Stop, Auto Reverse Run, Reciprocating



8 Memory Programs


Selecting each number brings up the stored torque, rotation speed, and rotation mode settings. You can save new settings in memory.



Continuous use


Available for 3 hours and 30 minutes with 1 hour and 30 minutes of full charge (※ Based on continuous no-load operation)






Reciprocating mode

The motor reciprocates in the forward and reverse directions.
Rotate in the normal forward 150 ° / reverse 30 °
the rotation value changes flexibly according to the load value.



Auto Reverse Stop mode


Set while the motor rotates in the forward direction reverse
when the torque value is reached rotate If the load is lower
than the set torque value the motor stops.



Auto Reverse Run mode

Set while the motor rotates in the forward direction reverse
when the torque value is reached rotate If the load is lower
than the set torque value the motor will rotate in the forward
direction again.





DRIVE HANDPIECE	
Model	EH-C500D
Battery type	Li-Polymer
Battery capacity	7.4V DC / 360 mAh
Speed range	120 ~ 500 rpm ±10%
Torque stages	0.5/1.0/1.5/2.0/2.5/3.0/3.5/4.0 N.cm ±10%
Dimension	27mm(W) x 23mm(L) x 157mm(H)
Weight	98g
CHARGING CRADLE	
Model	EH-CS12
Input	12V DC / 1.0A
Output	8.4V DC / 360 mAh
Charging time	Approx. 90min
Dimension	105mm(W) x 111mm(L) x 51mm(H)
Weight	80g
ADAPTOR	
Input	100 ~ 240 VAC, 0.3A
Output	12V DC / 0.5A
ANGLE	
Scerw driver	Shank Ø2.35mm, ISO1797-1 Type1
Chuck type	Push button chuck
Weight	50g

endoit[®]

Designed by **MICRO NX**



MICRO-NX Co., Ltd.

22, Maeyeo-ro 1-gil, 41059 Republic of Korea
Tel: 82 53 650 1000 Fax: 82 53 650 1001



MICRO-NX Co., Ltd. branch office

Karl-Marx-Str. 6 16540 Hohen Neuendorf Germany
Tel: 49 (0)3303 5412323 Fax: 49 (0)3303 5412324



Official Website



Official Youtube Chanel