



Please read this manual before operating

UDS-J ULTRASONIC SCALER INSTRUCTION MANUAL



www.glwoodpecker.com

GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

Contents

1. The installation and components of equipment	1
2. Installation and adjustment	2
3. Operation methods and function instruction	2
4. Cleaning, Disinfection and Sterilization	3
5. Contraindication	9
6. Storage and maintenance	9
7. Troubleshooting and notes	9
8. Precaution	13
9. Transportation	13
10. Working condition	13
11. After service	14
12. Symbol instruction	14
13. Environmental protection	15
14. Manufacturer's right	15
15. EMC-Declaration of conformity.....	15
16. Statement	18
TABLE OF OPERATING POWER OF THE TIPS.....	19

1. The installation and components of equipment

1.1 Instruction

Guilin Woodpecker Medical Instrument Co., Ltd. is a professional manufacturer in researching, developing and producing ultrasonic scalers. The product is mainly used for teeth cleaning and also an indispensable equipment for tooth disease prevention and treatment.

The ultrasonic scaler UDS-J has scaling, perio functions with the following features:

1.1.1 Automatic frequency tracking ensures that the machine always works on the best frequency and more steadily.

1.1.2 Digitally controlled, easy operation and more efficient for scaling.

1.2 Components

1.2.1 The components of machine are listed in the packing list.

The scaling tips and other accessories are not listed in this instruction manual completely. The detail can be found in the instruction for tips and packing list.

1.2.2 Product performance and structure

Ultrasonic scaler is composed of electro circuit, water way and ultrasonic transducer.

1.2.3 Scope of application

Ultrasonic scaler UDS-J is used for the dental calculus elimination.

1.3 The main technical specifications

1.3.1 Power supply input: 220-240V~ 50Hz/60Hz 150mA

1.3.2 Main unit input: 24V~ 50Hz/60Hz 1.3A

1.3.3 Output power: 3W to 20W

1.3.4 Output tip vibration frequency: 30kHz±3kHz

1.3.5 Output half-excursion force: <2N

1.3.6 Output primary tip vibration excursion: ≤100μm

1.3.7 Main unit fuse: T1.6AL 250V

1.3.8 Power supply fuse: T0.5AL 250V

1.3.9 Water pressure: 0.01MPa to 0.5MPa

1.3.10 Main unit weight: 0.8 kg

1.3.11 Power supply weight: 1 kg

1.3.12 Operating mode: Continuous operation

1.3.13 Type of protection against electric shock: class II equipment

1.3.14 Degree of protection against electric shock: Type BF applied part


1.3.15 Degree of protection against harmful ingress of water: Ordinary equipment (IPX0), the foot switch is drip-proof equipment (IPX1)

1.3.16 Applied part of the equipment: handpiece and tip

1.3.17 Degree of safety of application in the presence of a Flammable Anesthetic Mixture with air, Oxygen or Nitrous Oxide: Equipment not suitable for being used in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

1.4 Components instruction

The components of the equipment are as showed in picture 1.



Picture 1

2. Installation and adjustment

2.1 Open the packing box, make sure that all the parts and accessories are complete according to the packing list.

2.2 Take the main unit out of the box and put it on a stable plane, keep the main unit straight to the operator.

2.3 Turn the water control knob towards clockwise direction to the max and turn the potentiometer towards clockwise direction to a suitable position. [note2]

2.4 Insert the plug of the foot switch to its socket.

2.5 Connect one end of the water pipe to the water entrance, and the other end to the clean water source.

2.6 Connect the output end of power supply with main unit and get through to the power

2.7 Press the power switch of the main unit, then the power indicator shines.

3. Operation methods and function instruction

3.1 Direct the pit of the potentiometer at the "1" dial on the corer before turning on the scaler, make the main unit straight to the operator and turn the water control knob towards clockwise direction about three circles to the maximum.

3.2 The normal frequency is $30\text{kHz} \pm 3\text{kHz}$. With the high frequency, a light touch and a certain to-and-fro motion will eliminate the tartar without obvious heating. Overexertion and overstay are forbidden.

3.3 The way of assembly and disassembly of the scaling tips is as showed in picture 2.

3.4 The choice and operation methods of tips is shown in detail in attached materials of the equipment.

3.5 Vibrating intensity: Adjust the vibrating intensity

as you need, generally turn the knob to the

middle

grade. According to patient's different sensitivity and the rigidity of the gingival tartar, to adjust the vibrating intensity during the clinical treatment.

3.6 Water volume adjust: Step on the foot switch, and the tip begins to vibrate, then turn the water control knob to form spray, so as to cool down the handpiece and clean the teeth.

3.7 The handpiece can be handled in the same gesture as a pen in hand.

3.8 During the clinical treatment, be sure not to make the end of tip touch the teeth vertically

and not to make the tip overexert on the surface of the teeth, in case of hurting the teeth and damaging the tip.

3.9 After finishing operation, keep the machine working for 30 seconds with the water supply in order to clean the handpiece and the scaling tip.

3.10 Unscrew the scaling tip and sterilize it.



Picture 2:
Fasten the scaling tip by wrench with the hand between thumb and index finger.

4. Cleaning, Disinfection and Sterilization

The cleaning, disinfection and sterilization of tip, and wrench (include 1# torque wrench and Endo wrench) are as follow.

Unless otherwise stated, they will be hereinafter referred to as “products”.

Warnings

The use of strong detergent and disinfectant (alkaline pH>9 or acid pH <5) will reduce the life span of products. And in such cases, the manufacturer takes no responsibility.

This products shall not be exposed to high temperature above 138°C.

Processing limit

The products have been designed for a large number of sterilization cycles. The materials used in manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products. The maximum number of sterilizations for tips is 300 times. And for wrench, it is 1000 times.

4.1 Initial processing

4.1.1 Processing principles

It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic, especially with regard to the

additional requirements for the inactivation of prions.

4.1.2 Post-operative treatment

The post-operative treatment must be carried out immediately, no later than 30 minutes after the completion of the operation. The steps are as follows:

Tools: Endo wrench or 1# torque wrench, tray, clean and dry soft cloth

1. Let the UltrasonicScaler works for 20-30 seconds at maximum water volume to flush the handpiece and tip;

2. Soak the soft cloth with pure water (or distilled water or deionized water), and then wipe all the surfaces of the handpiece and tip until the surface of them is not stained;

3. Dry the handpiece and tip with a clean, soft cloth;

4. Remove the tip from handpiece with Endo wrench or 1# torque wrench provided by Guilin Woodpecker Medical Instrument Co., Ltd, and then put the tip and wrench into a clean tray.

Notes

a) The water used here must be pure water, distilled water or deionized water.

4.2Cleaning

The cleaning should be performed no later than 24 hours after the operation.

The cleaning can be divided into automated cleaning and manual cleaning. Automated cleaning is preferred if conditions permit.

4.2.1Automated cleaning

•The cleaner is proved to be valid by CE certification in accordance with EN ISO 15883.

▪There should be a flushing connector connected to the inner cavity of the product.

▪The cleaning procedure is suitable for the product, and the irrigating period is sufficient.

It is recommended to use a washer-disinfector in accordance with EN ISO 15883. For the specific procedure, please refer to the automated disinfection section in the next section "Disinfection".

Notes

a) The cleaning agent does not have to be pure water. It can be distilled water, deionized water or multi-enzyme. But please ensure that the selected cleaning agent is compatible with the product.

b) In washing stage, the water temperature should not exceed 45 °C, otherwise the protein will solidify and it would be difficult to remove.

c) After cleaning, the chemical residue should be less than 10mg / L.

4.3Disinfection

Disinfection must be performed no later than 2 hours after the cleaning phase. Automated disinfection is preferred if conditions permit.

4.3.1 Automated disinfection-Washer-disinfector

•The washer-disinfector is proved to be valid by CE certification in accordance with EN ISO 15883.

•Use high temperature disinfection function. The temperature does not exceed 134 ° C, and the disinfection under the temperature cannot exceed 20 minutes.

•The disinfection cycle is in accordance with the disinfection cycle in EN ISO 15883.

Cleaning and disinfecting steps by using Washer-disinfector

1. Carefully place the product into the disinfection basket. Fixation of product is needed only when the product is removable in the device. The products are not allowed to contact each other.

2. Use a suitable rinsing adaptor, and connect the internal water lines to the rinsing connection of the washer-disinfector.

3. Start the program.

4. After the program is finished, remove the product from the washer-disinfector, inspect (refer to section "Inspection and Maintenance") and packaging (refer to chapter "Packaging"). Dry the product repeatedly if necessary (refer to section "Drying").

Notes

a) Before use, you must carefully read the operating instructions provided by the equipment manufacturer to familiarize yourself with the disinfection process and precautions.

b) With this equipment, cleaning, disinfection and drying will be carried out together.

c) Cleaning: (c1) The cleaning procedure should be suitable for the product to be treated. The flushing period should be sufficient (5-10 minutes). Pre-wash for 3 minutes, wash for another 5 minutes, and rinse it for twice with each rinse lasting for 1 minute. (c2) In the washing stage, the water temperature should not exceed 45 °C, otherwise the protein will solidify and it is difficult to remove. (c3) The solution used can be pure water, distilled water, deionized water or multi-enzyme solution, etc., and only freshly prepared solutions can be used. (c4) During the use of cleaner, the concentration and time provided by manufacturer shall be obeyed. The used cleaner is neodisher MediZym (Dr. Weigert).

d) Disinfection: (d1) Direct use after disinfection: temperature ≥ 90 ° C, time ≥ 5 min or $A_0 \geq 3000$;

(d2) Sterilize it after disinfection and use: temperature ≥ 90 ° C, time ≥ 1 min or $A_0 \geq 600$

(d3) For the disinfection here, the temperature is 93 ° C, the time is 2.5 min, and $A_0 > 3000$

e) Only distilled or deionized water with a small amount of microorganisms (< 10 cfu/ml) can be used for all rinsing steps. (For example, pure water that is in accordance with the European Pharmacopoeia or the United States

Pharmacopoeia).

- f) After cleaning, the chemical residue should be less than 10mg / L.
- g) The air used for drying must be filtered by HEPA.
- h) Regularly repair and inspect the disinfectant.

4.4 Drying

If your cleaning and disinfection process does not have an automatic drying function, dry it after cleaning and disinfection.

Methods

1. Spread a clean white paper (white cloth) on the flat table, point the product against the white paper (white cloth), and then dry the product with filtered dry compressed air (maximum pressure 3 bar). Until no liquid is sprayed onto the white paper (white cloth), the product drying is completed.

2. It can also be dried directly in a medical drying cabinet (or oven). The recommended drying temperature is 80°C~120°C and the time should be 15~40 minutes.

Notes

- a) The drying of product must be performed in a clean place.
- b) The drying temperature should not exceed 138 °C;
- c) The equipment used should be inspected and maintained regularly.

4.5 Inspection and maintenance

In this chapter, we only check the appearance of the product.

4.5.1 Check the product. If there is still visible stain on the product after cleaning/disinfection, the entire cleaning/disinfection process must be repeated.

4.5.2 Check the product. If it is obviously damaged, smashed, detached, corroded or bent, it must be scrapped and not allowed to continue to be used.

4.5.3 Check the product. If the accessories are found to be damaged, please replace it before use. And the new accessories for replacement must be cleaned, disinfected and dried.

4.5.4 If the service time (number of times) of the product reaches the specified service life (number of times), please replace it in time.

4.6 Packaging

Install the disinfected and dried product and quickly package it in a medical sterilization bag (or special holder, sterile box).

Notes

- a) The package used conforms to ISO 11607;
- b) It can withstand high temperature of 138 °C and has sufficient steam permeability;
- c) The packaging environment and related tools must be cleaned regularly to ensure cleanliness and prevent the introduction of contaminants;
- d) Avoid contact with parts of different metals when packaging.

4.7 Sterilization

Use only the following steam sterilization procedures (fractional pre-vacuum procedure*) for sterilization, and other sterilization procedures are prohibited:

1. The steam sterilizer complies with EN13060 or is certified according to EN 285 to comply with EN ISO 17665;
2. The highest sterilization temperature is 138 °C;
3. The sterilization time is at least 4 minutes at a temperature of 132 °C / 134 °C and a pressure of 2.0 bar ~ 2.3 bars.
4. Allow a maximum sterilization time of 20 minutes at 134 °C.

Verification of the fundamental suitability of the products for effective steam sterilization was provided by a verified testing laboratory.

Notes

- a) Only products that have been effectively cleaned and disinfected are allowed to be sterilized;
- b) Before using the sterilizer for sterilization, read the Instruction Manual provided by the equipment manufacturer and follow the instructions;
- c) Do not use hot air sterilization and radiation sterilization as this may result in damage to the product;
- d) Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness.

*Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used here is to perform steam sterilization through three pre-vacuums.

4.8 Storage

4.8.1 Store in a clean, dry, ventilated, non-corrosive atmosphere with a relative humidity of 10% to 93%, an atmospheric pressure of 70KPa to 106KPa, and a temperature of -20 °C to +55 °C;

4.8.2 After sterilization, the product should be packaged in a medical sterilization bag or a clean sealing container, and stored in a special storage cabinet. The storage time should not exceed 7 days. If it is exceeded, it should be reprocessed before use.

Notes

- a) The storage environment should be clean and must be disinfected regularly;
- b) Product storage must be batched and marked and recorded.

4.9 Transportation

1. Prevent excessive shock and vibration during transportation, and handle with care;

2. It should not be mixed with dangerous goods during transportation.

3. Avoid exposure to sun or rain or snow during transportation.

The cleaning and disinfection of main unit and handpiece are as follows.

1 Pre-Op processing

Before each use, the handpiece and main unit must be cleaned and disinfected.

The specific steps are as follows:

1.1 Manual cleaning steps:

1.1.1 Wet the soft cloth completely with distilled water or deionized water, and then wipe all the surfaces of the handpiece and main unit until the surface of them is not stained.

1.1.2 Wipe the surface of the handpiece and main unit with a dry soft nap-free cloth.

1.1.3 Repeat the above steps at least 3 times.

Notes:

a) Use distilled water or deionized water for cleaning at room temperature.

1.2 Manual disinfection steps:

1.2.1 Soak the dry soft cloth with 75% alcohol

1.1.2 Wipe all the surfaces of the handpiece and main unit with a wet soft cloth for at least 3 minutes.

1.1.3 Wipe the surface of the handpiece and main unit with a dry soft nap-free cloth.

Notes:

a) The cleaning and disinfection must be performed within 10min before use.

b) The disinfectant used must be used immediately, no foaming is allowed.

c) In addition to 75% alcohol, you can use non-residue disinfectants such as Oxytech from Germany, but you must respect the concentration, temperature and time specified by the disinfectant manufacturer.

2 Post-Op processing

After each use, clean and disinfect the handpiece and main unit within 30 minutes. The specific steps are as follows:

Tools: Nap-free soft cloth, tray

2.1 Remove the tip from handpiece with Endo wrench or 1# torque wrench provided by Guilin Woodpecker Medical Instrument Co., Ltd, and then put the tip and wrench into a clean tray.

2.2 Soak the nap-free soft cloth with distilled water or deionized water, and then wipe all the surfaces of the handpiece and main unit until the surface of them is not stained.

2.3 Wet the dry soft cloth with 75% alcohol, and then wipe all surfaces of the handpiece and main unit for 3 minutes.

2.4 Put the handpiece and main unit back into the clean storage area.

Notes:

a) The cleaning and disinfection must be performed within 10min before use.

b) The disinfectant used must be used immediately, no foaming is allowed.

c) In addition to 75% alcohol, you can use non-residue disinfectants such as Oxytech from Germany, but you must respect the concentration, temperature and time specified by the disinfectant manufacturer.

5. Contraindication

5.1 The patient who has hemophilia is not allowed to use this equipment.

5.2 The patient or doctor who with heart pacemaker is forbidden to use this equipment.

5.3 The heart disease patient, pregnant woman and children should be cautious to use the equipment.

6. Storage and maintenance

6.1 The equipment should be handled carefully and lightly. Be sure that it is far from the vibration, and is installed or kept in a cool, dry and ventilated place.

6.2 Don't store the machine together with the articles that are combustible poisonous, caustic, or explosive.

6.3 This equipment should be stored in a room where the relative humidity is 10% ~ 93%, atmospheric pressure is 70kPa to 106kPa, and the temperature is -20°C ~ +55°C.

6.4 Please turn off the power switch and pull out the power plug when the equipment is not used. If the machine is not used for a long time, please make it get through to the power and water once per month for five minutes.

7. Troubleshooting and notes

7.1 Troubleshooting

Fault	Possible	Solutions
The scaling tip doesn't vibrate and there is no water flowing out when stepping on the foot switch.	The power line plug is in loose contact.	Make the plug insert to the socket well.
	The foot switch is in loose contact.	Insert the foot switch to its socket tightly.
	The fuse of transformer is broken.	Open the power box, change a new T0.5AL250V fuse.
	The fuse in the main unit is broken.	Take off the cover, change a new T1.6AL250V fuse.

Fault	Possible	Solutions
The scaling tip doesn't vibrate but there is water flowing out when stepping on the foot switch.	The tip hasn't been screwed on the handpiece tightly.	Screw the tip on the handpiece tightly (picture 2).
	The connect plug between the handpiece and the circuit board is in loose contact.	Contact our dealers or us.
	Something wrong with the handpiece	Contact our dealers or us [note 1].
The handpiece generates heat.	The water control knob is in a low grade.	Turn on the water control knob to a higher grade [note 2].
	The potentiometer is damaged [note 2].	Change a new one.
The scaling tip vibrates but there is no fine spray when stepping on the foot switch.	The water control knob is turn off.	Turn on the water control knob [note 2].
	There is impurity in the solenoid valve.	Clean inside of solenoid valve (picture 5).
	The water pipe is blocked.	Clean the water pipe by multi-function syringe [note 3].
There is still water flowing out after the power is off.	There is impurity in the solenoid valve.	Clean inside of solenoid valve (picture 5).
The amount of spouting water is too little.	The water control knob is in a low grade.	Turn on the water control knob to a higher grade [note 2].
	The water pressure is not high enough.	Make the water pressure higher.
	The water pipe is blocked	Clean the water pipe by multi-function syringe [note2].

Fault	Possible	Solutions
The vibration of the tip becomes weak.	The tip hasn't been screwed on to the handpiece tightly.	Screw the tip on the handpiece tightly (picture 2).
	The tip is vibrated loose.	Screw the tip on the handpiece tightly (picture 2).
	The tip is damaged [note3].	Change a new one.
The vibration is too strong and the potentiometer is failure.	The potentiometer is damaged [note 2].	Change a new one.

If the troubles still can't be solved, please contact with the local distributors or our company.

7.2 Notice

[**Note 1**] The disassembly of the handpiece (As showed in picture 3):

a) Remove the screw from the cover, pull out the potentiometer vertically, then take off the cover from the end of the machine lightly (there is a line connecting the cover and the machine, don't use too much strength). The inner of the machine is as showed in picture 3.

b) Pull out the water pipe in the handpiece cable from the coupling between the water control knob and the water pipe.

c) Pull out the lead plug from circuit board and untie the string.


d) Hold the joint of handpiece cable and the main unit and push it into the main unit about 1cm, then pull out the cutting sleeve from the cable.

e) Take off the handpiece from the main unit, and the disassembly is finished.

The assembly of handpiece is on the contrary. Be sure not to assemble the lead plug in a wrong direction, otherwise the tip will be electriferous.

Check-up method: get through to the power, step on the foot switch, and check-up the scaling tip with electric pen. If the electric pen shines, the scaling tip is electriferous. Turn off the machine and then insert the lead plug correctly.

[**Note 2**] Turn the water control knob towards anticlockwise direction, when the knob can't be turned any more, it comes to the min. On the contrary direction, the water volume increases step by step till the knob is back-out. The grade of the potentiometer is from gear 1 to 9. The ninth grade is the max. Be sure not to overdo.



Picture 3 The disassemblt of handpiece

[Note 3] To clean the water pipe with the multi-function syringe of the dental unit (as showed in the picture 4):

- Snip the water pipe at a distance of 10cm to 15cm from the water entrance.
- Turn on the power switch, get through to the power.
- Connect the multi-function syringe of the dental unit to the water pipe.
- Step on the foot switch.
- Turn on the switch of the multi-function syringe, press the air or water into the water pipe in the machine, then eliminate the impurity in the water pipe.



Picture 4 Multi-fuction syringe

[Note 4] If the scaling tip has been screwed on tightly and there is fine spray too, the following phenomena shows that the scaling tip was damaged:

- The vibrating intensity and the pulverization degree become weak obviously.
- When operating, there is some buzz when the scaling tip is working.

8. Precaution

- 8.1 Notice when using equipment
- 8.2 Keep the scaler clean before and after operation.
- 8.3 The scaling tip, wrench and handpiece must be sterilized before each treatment.
- 8.4 Don't screw the scaling tip when stepping on the foot switch.
- 8.5 The scaling tip must be fastened. There must be fine spray coming out from the tip when operating.
- 8.6 Change a new one when the tip is damaged or worn excessively.
- 8.7 While scaler working ,the heat of scaling tip will become higher if there is no water flowing out.Please keep the water flow smoothly.
- 8.8 Don't twist or rub the tip.
- 8.9 Don't use impure water source, and be sure not to use normal brine instead of pure water source.
- 8.10 If use the water source without hydraulic pressure, the water surface should be one meter higher than the head of the patient.
- 8.11 After operating, turn off electrical source, and then pull out the plug.
- 8.12 As a professional manufacturer of medical instruments, we are only responsible for the safety on the following conditions:
 - 8.12.1 The maintenance, repair and modification are made by the manufacturer or the authorized dealer.
 - 8.12.2 The changed components are original of "WOODPECKER" and operated correctly according to instruction manual.
- 8.13 Please put the power plug into the socket easy to pull out, to make sure it can be pull out in emergency.
- 8.14 The power supply is considered as a part of ME equipment.This device can only be equipped with the special power supply of Guilin Woodpecker Medical Instrument Co., Ltd.
- 8.15 The power supply is NOT waterproof. Please keep it dry and away from the water.
- 8.16 The screw thread of the scaling tips produced by other manufacturers maybe coarse, rusty and collapsed, which will damage the screw thread of the handpiece irretrievably. Please use "WOODPECKER" brand scaling tip.
- 8.17 Please select a suitable power when using different type of tips (refer to "TABLE OF OPERATING POWER OF THE TIPS").

9. Transportation

- 9.1 Excessive impact and shake should be prevented in transportation. Lay it carefully and lightly and don't invert it.
- 9.2 Don't put it together with dangerous goods during transportation.
- 9.3 Avoid solarization and getting wet in rain or snow during transportation.

10. Working condition


- 10.1 Environment temperature: +5°C to +40°C
- 10.2 Relative humidity: 30% ~75%
- 10.3 Atmosphere pressure: 70kPa to 106kPa
- 10.4 A temperature of the water at the inlet: not higher than +25°C

11. After service

We offer one year free repair to the equipment according to the warranty card.

The repair of the equipment should be carried out by professional technician. We are not responsible for any irretrievable damage caused by the not professional person.

Note: "P" was put on the valve seat to designate the water entrance.



Fasten the nut

Fasten the sunk screw

Water entrance

Picture 5 The assembly and disassembly of the solenoid valve

12. Symbol instruction



Trademark



Follow Instructions for Use

IPX1

Anti-drip device

IPX0

Ordinary equipment



Alternating current



Foot switch interface



Date of manufacture



Manufacturer



Class II equipment



Type BF applied part

~24V

~24V power supply input interface



Used indoor only



Recovery




Power switch



Appliance compliance WEEE directive





13. Environmental protection

Please dispose according to the local laws.

14. Manufacturer's right

We reserve the rights to change the design of the equipment, the technique, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some differences between blueprint and real equipment, take the real equipment as the norm.


15. EMC-Declaration of conformity

Guidance and manufacturer's declaration - electromagnetic emissions		
The model UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are intended for use in the electromagnetic environment specified below. The customer or the user of the model UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance & Declaration — electromagnetic immunity			
The models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are intended for use in the electromagnetic environment specified below. The customer or the user of the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for Input/output lines	±2kV for power supply lines ±1kV for interconnecting cable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED require continued operation during power mains interruptions, it is recommended that the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance & Declaration - Electromagnetic immunity

The models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are intended for use in the electromagnetic environment specified below. The customer or the user of the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3V 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 3V $d=1.2 \times P^{1/2}$ 80 MHz to 800 MHz $d=2.3 \times P$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz end 800 MHz. the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are used exceeds the applicable RF compliance level above, the model UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED			
The model UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED is intended for use in electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the models UDS-J, UDS-K, UDS-K LED, UDS-L, UDS-L LED, UDS-A, UDS-A LED, UDS-P, UDS-E, UDS-P LED, UDS-E LED, D1, D3, D5, D7, D3 LED, D5 LED, D7 LED are recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency m		
	150kHz to 80MHz $d=1.2 \times P^{1/2}$	80MHz to 800MHz $d=1.2 \times P^{1/2}$	800MHz to 2,5GHz $d=2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference Avoid using the device in high electromagnetic environment.

16. Statement

All rights of modifying the product are reserved to the manufacturer without further notice. The pictures are only for reference. The final interpretation rights belong to GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD. The industrial design, inner structure, etc, have claimed for several patents by WOODPECKER, any copy or fake product must take legal responsibilities.

TABLE OF OPERATING POWER OF THE TIPS

Scaling	
Tip Model	Power
G 1	1-9
G 2	1-9
G 3	1-9
G 4	1-9
G 5	1-9
G 6	1-9
G 7	1-9
G 8	1-9
G 9	1-9
G 10	1-9
G 11	1-9

Periodontics	
Tip Model	Power
P1	1-6
P2L	1-2
P2LD	1
P2R	1-2
P2RD	1
P3	1-3
P3D	1-3
P4	1-3
P4D	1-3

Endodontics	
Tip Model	Power
E1	-
E2	-
E3	-
E3D	-
E4	-
E4D	-
E5	-
E5D	-
E8	-
E9	-
E10	-
E10D	-
E11	-
E11D	-
E14	-
E15	-

Cavity Preparation	
Tip Model	Power
SB1	1-6
SB2	1-6
SB3	1-6
SBL	1-6
SBR	1-6

Scan and Login website
for more information



Guilin Woodpecker Medical Instrument Co., Ltd.
Information Industrial Park, Guilin National High-Tech
Zone, Guilin, Guangxi, 541004 P. R. China

Tel:

Europe Sales Dept.: +86-773-5873196

North/South America & Oceania Sales Dep.: +86-773-5873198

Asia & Africa Sales Dep.: +86-773-5855350 Fax: +86-773-5822450

E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com

Website: <http://www.glwoodpecker.com>