

# CuringPen-E

3 LEDs curing light

## **Technical Specifications:**

Dimension : 19 x 16.5 x 7.5 cm

Weight : 120 g
Lens Size : 8 mm
Power Supply :>1200

: >1200 times,

3000 mw/cm<sup>2</sup>, 3 Sec

Maximum Light Intensity : 3000 mw/cm<sup>2</sup> Spectral Range : 380-515 nm



#### 360° Rotatable Head

Easy to illuminate various treatment areas



#### Small Light Head

Easy to access to posterior molars which offers better patient experience



#### **Special Lens Design**

Unique lens design forms a concentrated beam for uniform curing



#### **Ultra Power**

- Light intensity now up to 3000mw/cm<sup>2</sup>
- High light intensity results in deeper cure depths which could curing resins more efficiently and with short time




### Broad Spectral Range 380 nm-515 nm

- Why broad spectral range?
- Most of resins use camphorquinone (CQ) materials, but some manufacturers choose alternative initiators such as TPO to reduce the yellow impact caused by CQ. TPO can only be activated by the light wave below 410nm. CuringPen-E has a broad spectrum which could effectively cure different resins.

#### **Working Modes**

•P1



Normal Mode

Light Intensity : 1200mw/cm<sup>2</sup>

Time Setting : 5s, 10s, 15s, 20s, 25s, 30s, 35s, 40s

•P2

High Mode

□ PB ② □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □
 □

•P3

Light Intensity : 2000mw/cm² Light Intensity : 3000mw/cm² Time Setting : 1s, 3s, 5s, 10s Time Setting : 1s, 3s, 5s

•P4

Eighteeth



Pulse Mode

Light Intensity : 1200mw/cm²
Time Setting : 5s, 10s, 15s, 20s


•P5




Ramp Mode

Light Intensity : 1200mw/cm²
Time Setting : 5s, 10s, 15s, 20s

Pulse mode



Ramp mode



**Eighteeth**